



Gel Distal End Pad

Worn inside a gel liner, these distal end pads are designed to provide added distal cushioning and support. Made of pure, soft silicone gel.

Order No.	Diameter
ESP-DPXS	4.45cm
ESP-DPS	5.08cm
ESP-DPM	5.39cm
ESP-DPL	5.87cm
ESP-DPXL	6.98cm

Silicone Gel Disc

Used inside a gel liner for added cushioning, compression or relief from irritation due to irregular limb shape or volumetric changes. 4mm thick with tapered edges.

Order No.	Diameter
ESP-GD	4"

Silicone Gel Strips

Used inside a gel liner for added cushioning, compression or relief from irritation due to irregular limb shape or volumetric changes. 2mm or 4mm thickness.

Order No.	W	L
2mm		
ESP-2GSA	1"	3"
ESP-2GSB	1"	4.5"
ESP-2GSC	1.25"	3"
ESP-2GSD	1.25"	4.5"
4mm		
ESP-4GSA	1"	3"
ESP-4GSB	1"	4.5"
ESP-4GSC	1.25"	3"
ESP-4GSD	1.25"	4.5"