




RUSH Foot® Product Part Number Guide



Product Code

Chopart = **CHO**
 HiPro = **HiP**
 LoPro = **LoP**
 RAMPAGE = **RAM**
 HiPro H2O = **H2H**
 LoPro H2O = **H2L**
 ROGUE H2O = **H2R**
 ROGUE = **ROG**
 * Foot Shell = **FS**
 ** Spectra Sock = **SS**
 EVAQ8 ROGUE = **EVQR**
 EVAQ8 HiPro = **EVQH**
 EVAQ8 LoPro = **EVQL**
 EVAQ8 = **EVRB**
 Rebuild Kid
 *** EVAQ8 Release Valve = **EVRV**

Activity Level

Low = **L**
 Moderate = **M**
 High = **H**
 High Impact = **P**

Foot Size

22 cm = 22
 23 cm = 23
 24 cm = 24
 25 cm = 25
 26 cm = 26
 27 cm = 27
 28 cm = 28
 29 cm = 29

Proximal Adapter

Static Male – LoPro®, HiPro®, HiPro H2O, RAMPAGE, ROGUE®, ROGUE H2O, EVAQ8 = **MS**
 Static Female - HiPro®, HiPro H2O only, and EVAQ8 = **FS**
 Rotatable Male - HiPro®, HiPro H2O only, and EVAQ8 = **MR**
 Rotatable Female - HiPro®, HiPro H2O only, and EVAQ8 = **FR**

*If not selected, the proximal adapter will default to "MS", Static Male

CATEGORY CHART

Weight lbs	90-116	117-146	147-180	181-218	219-261	262-310	311-360
Low Activity	1	2	3	4	5	6	7
Mod Activity	1	2	3	4	5	6	7
High Activity	2	3	4	5	6	7	8
High Impact	2-3	3-4	4-5	5-6	6-7	7-8	8-9

*If not selected, the stiffness will default based on weight and activity.

Patient Weight

90-99 lbs = **09**
 100-109 lbs = **10**
 110-119 lbs = **11**
 120-129 lbs = **12**
 130-139 lbs = **13**
 140-149 lbs = **14**
 150-159 lbs = **15**
 160-169 lbs = **16**
 170-179 lbs = **17**
 180-189 lbs = **18**
 190-199 lbs = **19**
 ...continue pattern
 290-299 = **29**
 ...continue pattern
 390-399 = **39**

Activity Level Definitions

Low - Golfing, Gardening, Walking
 Moderate - Jogging, Hiking, Skiing
 High - Baseball, Running, Wakeboarding
 High Impact - Sprinting, Basketball, Dirt Bike Riding


Foot Shell

Left Light = **L1**
 Left Dark = **L3**
 Right Light = **R1**
 Right Dark = **R3**

* **Customers ordering foot shell only will indicate foot shell number only**
 ** **Customers ordering spectra sock only will indicate spectra sock number only**
 *** **RUSH EVAQ8 comes with a standard installation kit. Valve not included, but available for additional cost.**

PART NUMBER EXAMPLE

An amputee weighing **225lbs** wears a **28 cm** foot and is considered a **high impact user**. The prosthetist would like to purchase a HiPro® foot with **Static Male proximal adapter** and a **Right Dark foot shell**.



Cascade Orthopedic Supply, Inc.
 Toll Free: 800-888-0865
 Fax: 800-847-9180
 www.cascade-usa.com

